TG Cloud App Manual

General Description

The app (Android only) is designed to work with our Data Capture 7 software. Data Capture 7 running on one of your network computers can monitor your Temperature Guard unit and send data to a server in the cloud. The app can then retrieve the last ten readings from the cloud server and graph them.

The name of the unit is determined by Data Capture 7. Once you check the checkbox and click Save, Data Capture will automatically create a new file on the cloud website and retrieve a unique 32 digit code.

You will need to enter that code into the app. We recommend you email it to yourself as typing it in manually may result in errors.

Add/Modify Units	
Save Add New Edit IP Set	Data Logs Print Help
 Temperature Guard Units CM705E_R25 @ 10.0.0.76 CVG 41-2 @ 10.0.0.80 M308R10 @ 10.0.0.210 M705E_v128 @ 10.0.0.199 Rev 20 @ 10.0.0.188 VM505E @ 10.0.0.109 VM608E_TEST_UNIT @ 10.0.0.183 	Uni: 9 IP Address 10.0.0.210 Test Name of Device M308: Model M308: netGuard Send alarm emails to: Email addresses for this unit. One email per line. I Collect Data from this Unit Temperature Measuring Units (C or F) C Allow the data collected by this unit to be available to 1 G Cloud App Host ID 19.73b1a9ebaaba224086ad367 Date/Time of last communication 4/29/2025 3.05:23 PM Enter notes to include with alarm emails Enter notes to include with alarm emails

Adding a Unit

When the app is first run you will need to add at least one unit. Select the Add Unit menu item from the Android menu (three dots).

Paste in the 32 digit code generated by Data Capture.

Tap the Save button.



The app will communicate with the website and save the unit into the app.

If the name comes up as "dUN" do not worry it will automatically update when there is data on the website.



To delete a unit just select it from the drop down menu and select Delete Unit.

Select OK to delete, Cancel to not delete.

Using the app

To get data from your unit just select the name of the unit from the drop down list.

The app will communicate with the website and fill the chart with the latest data. Date and time is shown above the chart. Tap any sensor and the graph will appear. Graphs are best viewed in landscape mode.

10:55 🖱 🕯 🔊	🕿 🖌 🕯 100%	11:12 O 🔊		₹⊿	0 99%	10:39 🖱 🔊		₹⊿	97%
Home	:	Home			÷	Home			
Available Units M308R10	×	Available Units M705E_v128		•		Available Units M705E_v128		•	
M308R10									
M705E_v128									
	00:00 AM		05/02/2025	11:10 AM			05/02/2025	10:35 AM	
° Sensor	Data Status TOL	° Sensor	° Data	Status	TOL	° Sensor	° Data	Status	TOL
		Freezer VFC	2.9	ок	0	Freezer VFC	0.7	ок	0
		Refrigerator	36.6	ок	0	Refrigerator	39.1	ок	0
		Room Temp	69.9	ок	0	Room Temp	69.7	ок	0
		1K Resistor	30.6	ок	0	1K Resistor	30.6	ок	0
		Room Temperature	75.2	ок	0	Room Temperature	75.3	ок	0
		Room Humidity	39	ОК	0	Room Humidity	38	ок	0
		Freezer Door	CLOSED	Normal	0	Freezer Door	CLOSED	Normal	0
		Fridge Door	CLOSED	Normal	0	Fridge Door	CLOSED	Normal	0
		Backup Battery	3.9	OK	0	Backup Battery	3.9	ок	0
						44 43 42 41	Refrigerator		
•	• •		•			•	•		
	10:39 🔿 🔊					₹ ⊿ 197	7%		



If you have door, water, or other dry contact sensors you can tap and graph those as well. A second chart will replace the sensor chart. Tap in the chart to get back to the previous sensors chart. Scrolling is allowed in any of the white space on the screen.

			≇ ⊿ 🛙 100%
Home			÷
vailable Unit	5		
M705E_v1	28		•
	05	5/02/2025 11:35	АМ
reezer Door			
Date and	Time	Data	
5/02/2025 11:	35 AM	CLOSED	
5/02/2025 11:	30 AM	CLOSED	
5/02/2025 11:	30 AM	CLOSED	
5/02/2025 11:	30 AM		
,,	So Am	CLOSED	
5/02/2025 11:	30 AM	CLOSED	
5/02/2025 11: 5/02/2025 11:	30 AM 30 AM	CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM 70 AM	CLOSED CLOSED CLOSED CLOSED CLOSED CLOSED	
5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11: 5/02/2025 11:	30 AM 30 AM 30 AM 30 AM 30 AM 30 AM	CLOSED CLOSED CLOSED CLOSED CLOSED CLOSED	

Any issues please contact support@temperatureguard.com