



Quotation

Quote Number: Q1881

by Microtechnologies, Inc.

44 Washington Avenue, Suite 8
Berlin, CT 06037
1-888-248-7103 (Inside USA)
(860) 829-2710 (Outside USA)
(860) 356-4823 (Fax)
www.TemperatureGuard.com
sales@TemperatureGuard.com

Requested By:

Springfield Hospital
John Smith
123 Main Street
Berlin CT 06037

Phone
Fax
Email

Date: 5/31/2007 Quote valid for thirty days. Customer Reference Number:

ITEMS

Qty	Part No.	Description	Price Each	Line Total
1	MOLLY303	Ethernet Connectivity Refrigerator/Freezer Monitor	\$297.95	\$297.95
2	ACI/1K-2W-BP-V	(-55°F to 96°F) Temperature Sensor Probe in Vial	\$29.95	\$59.90
2	DR-1035W	Adhesive-Backed Plastic Door Sensor	\$19.95	\$39.90
1	SOFT-DC	Data Capture Desktop Software (Free Download)	\$0.00	\$0.00

Est. Shipping Weight: 3.5 lbs. **Shipping Method:** UPS Ground \$12.50 **\$410.25**

Terms: PO or Major Credit Card

ADDITIONAL COMMENTS

This quote is a complete system that allows Springfield Hospital to monitor the temperature inside both compartments (refrigerator and freezer) of a full size refrigerator.

The temperature sensors will go inside the refrigerator and freezer compartments and relate the temperature to the MOLLY303.

The door sensors will go on the doors of the refrigerator and freezer compartment and relate their open/closed state to the MOLLY303.

The Data Capture software will run on a computer that is on the network and collect the temperature readings and door status from the MOLLY303.

If either temperature goes out of range for too long or the doors are left open for too long, alarm email messages go out to a list of recipients and the MOLLY303 will alarm locally with its alarm buzzer.